

Data Sheet

QT PEROXY PURE



6% Hydrogen Peroxide Solution

QT Peroxy-Pure is a 6% Hydrogen Peroxide solution manufactured using USP water for injection (WFI). QT Peroxy-Pure is specifically formulated for removal of residues and cleaning of hard non-porous surfaces in clean rooms or other critical environments including, floors, walls, equipment, carts and tables. Bottles are filled and sealed aseptically in an ISO class 5 environment.

Each bottle is lot number and expiry date coded for complete traceability and straightforward record keeping. Bottles come fully assembled with a non-venting, ergonomic trigger sprayer with mist. The closed system (Bag in a Bottle) spray assembly maintains sterility throughout the products use. All sterile product is individually double-bagged in heat sealed poly bags and packaged in cases with individually sealed double case liners.

Features

- Contains 6% Hydrogen Peroxide with WFI
- 0.22 micron filtered
- Aseptically filled and sealed
- Biodegradable – completely decomposing into water and oxygen
- No VOCs – reduces VOC regulatory and compliance concerns
- Double-bagged and double sealed polybag case liner
- Lot specific records for each shipment:
 - Certificate of Analysis
 - Indelible lot & expiry markings on each container

Benefits

- Removes surface contaminants
- Low residue upon evaporating/no rinsing
- Pre-attached and ready to use trigger sprayer
- Convenient ready-to-use solutions require no mixing or activating

Applications

- Surface cleaning, walls, floors, tables
- Equipment cleaning



* Non-Sterile product will **NOT** include Gamma Indicators or Double-bagged.

| Part No. | Description | Quantity | Dimension | Weight | HS Code | Shelf Life |
|----------------|--|-----------------|------------------------------------|----------|---------------|------------|
| QT4500-16 | 6% Hydrogen Peroxide Bag in Bottle – 473mL – Sterile | 12 x 500mL/case | 350mm (L) x 280 mm (W) x 248mm (H) | 7.0 kg | 2847.00.10 00 | 24 months |
| QT4510-32 | 6% Hydrogen Peroxide Bag in Bottle – 946mL – Sterile | 6 x 1L/case | 305mm (L) x 215mm (W) x 320mm (H) | 8.0 kg | 2847.00.10 00 | 24 months |
| QT4501-16 | 6% Hydrogen Peroxide Trigger Spray – 473mL - Sterile | 12 x 500mL/case | 350mm (L) x 280mm (W) x 248mm (H) | 7.0 kg | 2847.00.10 00 | 24 months |
| QT4500-16NS | 6% Hydrogen Peroxide Trigger Spray – 473mL | 12 x 500mL/case | 350mm (L) x 280mm (W) x 248mm (H) | 7.0 kg | 2847.00.10 00 | 24 months |
| QT4511-32 | 6% Hydrogen Peroxide Trigger Spray – 946mL - Sterile | 6 x 1L/case | 305mm (L) x 215mm (W) x 320mm (H) | 8.0 kg | 2847.00.10 00 | 24 months |
| QT4510-32NS | 6% Hydrogen Peroxide Trigger Spray – 946mL | 6 x 1L/case | 305mm (L) x 215mm (W) x 320mm (H) | 8.0 kg | 2847.00.10 00 | 24 months |
| QT4520-128 | 6% Hydrogen Peroxide Screw Cap – 1 gal – Sterile | 4 x 1 gal/case | 410mm (L) x 260mm (W) x 300mm (H) | 16.26 kg | 2847.00.10 00 | 24 months |
| QT4520-128NS | 6% Hydrogen Peroxide Screw Cap – 1 gal | 4 x 1 gal/case | 410mm (L) x 260mm (W) x 300mm (H) | 16.26 kg | 2847.00.10 00 | 24 months |
| QT4530-9x20 | Presat Polyester Wipes 6% Hydrogen Peroxide, 9"x9" – Sterile | 10 pouches/case | 410mm (L) x 460mm (W) x 300mm (H) | 5.5 kg | 2847.00.10 00 | 12 months |
| QT4530-9x20NS | Presat Polyester Wipes 6% Hydrogen Peroxide, 9"x9" | 10 pouches/case | 410mm (L) x 460mm (W) x 300mm (H) | 5.5 kg | 2847.00.10 00 | 12 months |
| QT4530-12x20 | Presat Polyester Wipes 6% Hydrogen Peroxide, 12"x12" – Sterile | 10 pouches/case | 410mm (L) x 460mm (W) x 300mm (H) | 5.5 kg | 2847.00.10 00 | 12 months |
| QT4530-12x20NS | Presat Polyester Wipes 6% Hydrogen Peroxide, 12"x12" | 10 pouches/case | 410mm (L) x 460mm (W) x 300mm (H) | 5.5 kg | 2847.00.10 00 | 12 months |

PRODUCTS AVAILABLE STERILE OR NON-STERILE.