

## 3514500 Clino® CR1 EM-GMP Cleaning trolley for Cleanroom

Individualization  
on request!



### Product description

Cleaning trolley with EasyMop GMP® pre-preparation method and disposal for highly sensitive areas made of stainless steel.

- o Maximum efficiency and process reliability in a single-stage process
- o Made from the highest quality stainless steel
- o High-end surface quality through electropolishing
- o Touch-free operating concept



|     |   |   |     |   |   |   |
|-----|---|---|-----|---|---|---|
| ISO | 3 | 4 | 5   | 6 | 7 | 8 |
| GMP |   |   | A/B |   | C | D |

### Product specifications

|                                       |  |
|---------------------------------------|--|
| Article name   Article number         | Clino® CR1 EM-GMP   3514500  |
| Weight                                | 15.6 kg  |
| Working dimension                     | 61 x 53 x 70 cm  |
| Resting dimension                     | 61 x 53 x 70 cm  |
| Basic components                      | Basic module (stainless steel suitable for cleanroom use)<br>High-end finish<br>Autoclavable special rollers 100 mm<br>EasyMop GMP® system boxes with Mop Jump System<br>Device holder |
| Work area                             | Unit 1 x 35 L  |
| Disposal area                         | The cleaning trolley should be completely cleaned after use and stored in a dry place.   |
| Care instructions                     | The trolley's autoclave suitability depends on its individual components.  |
| Sterility                             |  |
| Pieces per PackingUnit   Selling unit | 1 Pieces   1 PU  |
| Pieces per box                        | 1 Pieces   |
| Cartons per pallet   Pallet           | 2   80 x 123 x 182 cm  |
| Customs tariff number                 | 87168000   |
| Country of origin                     | Germany  |
| EAN-Number                            | 4065237109007  |
| Compatible accessories                | Our sales department will be happy to advise you on suitable accessories   |

### Industry Sectors

### Assembly Hygiene



This classification is a non-binding recommendation. The user must check the individual suitability and compatibility for his individual area of application. The condition of the material must be checked regularly; damaged products must be replaced immediately. All information without guarantee. Subject to changes and deviations.



## 3514500 Clino® CR1 EM-GMP Cleaning trolley for Cleanroom

|    |         |                                      |                                      |
|----|---------|--------------------------------------|--------------------------------------|
| 1x | 3530171 | Module CR 1/2 - Facelift             | Autoclavable (121 °C for 20 minutes) |
| 1x | 3530142 | Bag support for CR1/CR4              | Autoclavable (121 °C for 20 minutes) |
| 4x | 3500100 | Stainless steel swivel castor 100 mm | Autoclavable (121 °C for 20 minutes) |
| 1x | 3530154 | Adapter for CR4                      | Autoclavable (121 °C for 20 minutes) |
| 2x | 3500200 | EasyMop® GMP system box              | Autoclavable (121 °C for 20 minutes) |
| 2x | 3500210 | MopJump®                             | Autoclavable (121 °C for 20 minutes) |
| 1x | 2017031 | Device holder 22 mm                  |                                      |
| 1x | 0899100 | Cardboard small                      |                                      |

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany  
www.pfennig-reinigungstechnik.com | info@pfennig-reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 16.03.2026;

